
TATA KELOLA SAMPEL LABORATORIUM PASIEN COVID-19

Dr. Hessyani P.T. Raranta, Sp.PK

CORONA VIRUS

Outbreak (Human)

MERS-CoV

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Tracheal aspirates are valuable specimens

- A reliable tracheal aspirate Gram stain can be used to direct initial empiric antimicrobial therapy and may increase the diagnostic value of the CPIS
- A negative tracheal aspirate (absence of bacteria or inflammatory cells) in a patient without a recent (within 72 hours) change in antibiotics has a strong negative predictive value (94%) for VAP and should lead to a search for alternative sources of fever


A circular inset illustration showing a patient's head and neck in profile. A bronchoscope is inserted into the mouth and down the trachea. A small tube is also shown inserted into the trachea, likely for aspirating secretions. The trachea and bronchi are highlighted in blue and red.

Table 1. Specimens to be collected from symptomatic patients

Guidance on specimen collection

| Specimen type | Collection materials | Transport to laboratory | Storage till testing | Comment |
|---|--|-------------------------|--------------------------------------|--|
| Nasopharyngeal and oropharyngeal swab | Dacron or polyester flocked swabs* | 4 °C | ≤5 days: 4 °C >5 days: -70 °C | The nasopharyngeal and oropharyngeal swabs should be placed in the same tube to increase the viral load. |
| Bronchoalveolar lavage | sterile container * | 4 °C | ≤48 hours: 4 °C >48 hours: -70 °C | There may be some dilution of pathogen, but still a worthwhile specimen |
| (Endo)tracheal aspirate, nasopharyngeal aspirate or nasal wash | sterile container * | 4 °C | ≤48 hours: 4 °C >48 hours: -70 °C | |
| Sputum | sterile container | 4 °C | ≤48 hours: 4 °C >48 hours: -70 °C | Ensure the material is from the lower respiratory tract |
| Tissue from biopsy or autopsy including from lung | sterile container with saline | 4 °C | ≤24 hours: 4 °C >24 hours: -70 °C | |
| Serum (2 samples acute and convalescent possibly 2-4 weeks after acute phase) | Serum separator tubes (adults: collect 3-5 ml whole blood) | 4 °C | ≤5 days: 4 °C >5 days: -70 °C | Collect paired samples: • acute – first week of illness • convalescent – 2 to 3 weeks later |
| Whole blood | collection tube | 4 °C | ≤5 days: 4 °C >5 days: -70 °C | For antigen detection particularly in the first week of illness |
| Urine | urine collection container | 4 °C | ≤5 days: 4 °C >5 days: -70 °C | |

Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases
Interim guidance
14 January 2020

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


Bio bottle :
For carry biological samples and specimens in a safe way for all transport modes,

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